

## Aidacare FL250 Bed End Set



- The FL250 Bed incorporates a sophisticated range of finishing options
- From our contemporary stock selection and made to order styles through to completely bespoke, any interior design desires can be met
- Attendant Control Panel - Allows the caregiver to lock out certain functions and also ensures bed operation is not just restricted to the handset

Made to order options include:

- **Tallo** - Timber or laminate surround and recessed timber veneer inner with a straight gap detail
- **Eyre** - Timber or laminate surround and recessed timber veneer inner with a curved gap detail.
- **Como** - Timber or laminatesurround and recessed timber veneer inner.
- **Sophia** - Timber surround with upholstered inner detail on headboard and recessed timber veneer detail on footboard (please contact us for custom orders)

- ✔ Contracted DVA product
- ✔ Contract approved



Code	Display Name
BEB046200	Bed End Set - FL250 - Single - Tahoe (Sublime Teak) - includes: Ctrl Panel + Hardware (stock)
BEB046210	Bed End Set - FL250 - Single - Vail (Milano Walnut) - includes: Ctrl Panel + Hardware (stock)
BEB046220	Bed End Set - FL250 - King Single - Tahoe (Sublime Teak) - includes: Ctrl Panel + Hardware (stock)
BEB046230	Bed End Set - FL250 - King Single - Vail (Milano Walnut) - includes: Ctrl Panel + Hardware (stock)
BEB046406	Bed End Set - FL250 Single - Como Straight with Radius Corners - Laminate - with NCP Cut-out - Sublime Teak

Code	Display Name
BEB046407	Bed End Set - FL250 King Single - Como Straight with Radius Corners - Laminate - with NCP Cut-out - Sublime Teak
BEB046408	Bed End Set - FL250 Single - Como Straight with Radius Corners, Select Beech includes NCP + Hardware
BEB046409	Bed End Set - FL250 King Single - Como Straight with Radius Corners, Select Beech includes NCP + Hardware