

Ceiling Hoist - Liko - Multirall 200



- The Multirall 200 hoist is a portable ceiling hoist that can also be used for room-to-room patient transfers
- Can be used to transfer a patient between rooms in a sitting position, easily moving from one rail system to another with the use of a transfer strap
- Its portability makes it an excellent and economical solution for nursing home and hospital care situations where lifting needs are temporary
- Can be easily attached or removed with no tools and minimum effort from ceiling rails to facilitate convenient transport between rooms
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Can be mounted two ways for ease of use

Warranty Information

- Frame: 10 Years
- Electrical Parts: 2 Years

- ✓ Contracted DVA product
- ✓ Contract approved

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Height	Safe Working Load
LIC359100	Multirall 200 - Portable Ceiling Hoist	Up to 1600 mm	200 kg