

ADVANCED ELECTRIC STRETCHER

Aspire for... Response, Support & Transfers



FEATURES & BENEFITS

- **X-ray capability and electric operation:** Radiolucent platform with electronic adjustment facilitates full body x-ray imaging.
- **Infection control:** Wipe clean surfaces adhere to stringent infection control procedures.
- **Medical Grade Castors:** 200mm Tente castors with central locking levers for optimal durability, traction and safety.
- **Grip and Glide:** Ergonomic grip handles and fifth centre wheel, enhance maneuverability and steering in restricted spaces and corridors.
- **Electrically adjustable positioning:** Improves patient independence and ease of use for staff.
- **Fold down side rails:** For safety and ease of transfers.
- **Comfortable pressure reducing mattress:** 100mm mattress helps to prevent shear and promote optimal skin microclimate.
- **Safety straps:** Integrated patient safety straps for peace of mind.
- **Standard inclusions:** Removable oxygen tank holder and IV pole.

**2 YEARS
WARRANTY**



CODE	DESCRIPTION	OVERALL LENGTH	OVERALL WIDTH	OVERALL HEIGHT	PATIENT SURFACE LENGTH	PATIENT SURFACE WIDTH	BACKREST ANGLE	TRENDELENBURG	CASTOR SIZE
TSE674230	Aspire Advanced Electric Stretcher	2140 mm	840 mm	600 - 910 mm	1920 mm	600 mm	0° to 70°	-12° to 12°	200mm

www.aspirecare.com.au
1300 133 120

