





BED MEASUREMENTS GUIDE




Relax Bed



Care Bed



Activ Care Bed



**Long single width, please see over the page for king single, long single and queen widths.*

BED MEASUREMENTS GUIDE

The Aspire ComfiMotion collection offers a range of beds and mattress in a variety of sizes. The dimensions listed below are for quick reference.

BED	LENGTH	WIDTH	HEIGHT RANGE		
			Activ Care	Care	Relax
Long Single	2030mm	920mm	195 – 700mm	280 – 650mm	280mm
King Single	2030mm	1070mm	195 – 700mm	280 – 650mm	280mm
Long Double	2030mm	1370mm	195 – 700mm		280mm
Queen	2030mm	1520mm	195 – 700mm		280mm

MATTRESS	LENGTH	WIDTH	HEIGHT
Long Single	2030mm	920mm	200 – 240mm
King Single	2030mm	1070mm	200 – 240mm
Long Double	2030mm	1370mm	200 – 240mm
Queen	2030mm	1520mm	200 – 240mm



comfimotion.com.au