

FURNITURE *Catalogue*



ATAMA  LIVING
by aidacare

ATAMA  LIVING
by aidacare

Atama has been a provider of customised furniture solutions to commercial and healthcare customers for over 30 years. We strive to balance quality and aesthetics with the functional outcomes required in modern healthcare environments.

Our recent acquisition by Aidacare brings an exciting new chapter as we expand and uplift our range and our ability to meet our customers needs across Australia.

Design, Manufacture, Deliver

While you explore this catalogue, rest assured all items can be tailored to suit. We will consult with you and strive to understand the details required to achieve the designs you desire. Whether we manufacture to suit or source and customise as required, we always ensure the highest quality standards are met. Your new furniture will be protected and delivered on time with the utmost care.

The Projects by Aidacare Display Centre also enables you to experience first hand all of the modern healthcare furnishings and equipment for your Greenfield and Brownfield developments and commercial projects.



FRAMES & UPHOLSTERY

Custom Frame Options _____ 05
 Custom Upholstery Options _____ 06

SEATING

Lounge & Highback Armchairs _____ 09
 Armchairs _____ 18
 Side Chairs _____ 25
 Day Chairs _____ 26
 Recliners _____ 28

TABLES

Coffee Tables _____ 33
 Side Tables _____ 35
 Dining Tables _____ 37

STORAGE

Sideboards & Buffets _____ 41
 Cabinets _____ 42
 • Bedside Cabinets
 • Storage Cabinets
 • TV and Media Cabinets

BEDS

Floorline - Aidacare FL250 _____ 46
 Adjustable - Aspire ComfiMotion _____ 48

ALFRESCO FURNITURE

_____ 50

WARRANTY

_____ 51

PROJECTS & FITOUTS

Projects by Aidacare Display Centre _____ 52

Frames & Upholstery

- CUSTOM FRAME OPTIONS
- CUSTOM UPHOLSTERY OPTIONS

CUSTOM FRAME OPTIONS

Metal FRAME

STAINLESS STEEL
BRASS
COPPER
ALUMINIUM

FINISHES



Wood GRAIN*

BEECH
VICTORIAN ASH
TASMANIAN OAK
EUROPEAN OAK
AMERICAN OAK

PREMIUM STAINS & LAMINATE



SELECT STAINS & LAMINATE



* All Select and Premium stains can be applied to any of the wood grains listed.

CUSTOM UPHOLSTERY OPTIONS

WARWICK ANDREAS



ACACIA BIRCH BRICK CHARCOAL CHOCOLATE PUMICE




DRIFTWOOD GLACIER INDIGO MIST OLD-ROSE REEF

WARWICK BRONX



ASH CHARTREUSE DUCK-EGG FROST JUNGLE LINEN STONE



MINK NAVY PAPRIKA PEPPER SLATE STEEL TURQUOISE

WARWICK OTWAY



ALMOND ASH CHARCOAL GRASS NAVY

MATERIALISED PRINTS



ACACIA RAIN GRAIN PRISCILLA BOUQUET NAVY ACORN WASABI SIENNA AMBER TANGIER FOG JAVA OLIVE MALUA MALLEE

WARWICK LUSTRELL CHARISMA VINYL



AMETHYST BLACK BUFFALO CABERNET CAMEL CHILI CLOUD



COBALT DOESKIN DOVE EARTH ECRU ELEPHANT FROST




GRASS GUNMETAL KIWI LEMON LIME NAVY OCEAN




OLIVE ORANGE PACIFIC PARCHMENT PECAN POLAR POPPY




PUMICE SADDLE SCARLET SEAFOAM SILVER SKY SLATE



SMOKE SNOW SUNSHINE TAUPE TEAL



STAINED GLASS APRICOT MILLA ORANGE



Seating

- LOUNGE & HIGHBACK ARMCHAIRS
- ARMCHAIRS
- SIDE CHAIRS
- DAY CHAIRS
- RECLINERS



Lounge & Highback Armchairs

All Lounge and Highback Armchairs can be manufactured in customised upholstery and frame finishes. See pages 4 to 7.

Images shown are for illustrative purposes only and may vary from custom options available.

AVAILABLE CONFIGURATIONS

KEY



Rosebud
Lounge Chair



Surrey
Lounge Chair



Sapphire
Lounge Chair



Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100350	Rosebud	720	740	955	540	520	490	520
FUL100290	Surrey	700	840	1050	500	500	490	600
FUL100330	Sapphire	700	740	880	500	490	490	400



Kew Mid Back
Lounge Chair



Elm
Lounge Chair



Banksia
Lounge Chair



Kew
Lounge Chair



Jordan
Lounge Chair



London
Lounge Chair



Cornwall
Lounge Chair



Gorden
Lounge Chair



Leeton
Lounge Chair



Bristol
Lounge Chair



Somerset
Lounge Chair



Devon
Lounge Chair



Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100180	Kew Mid Back	650	840	940	500	490	490	530
FUL100170	Elm	700	840	1000	510	470	490	550
FUL100120	Banksia	700	840	1000	510	400	490	550
FUL100100	Cornwall**	790	840	1040	480	480	490	600
FUL100240	Gorden*	700	840	980	510	490	490	530
FUL100270	Leeton	720	780	800	530	520	490	370

*Project only item: minimum order quantity required. **Alternative leg options available upon request.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100200	Kew	650	840	1040	500	490	490	600
FUL100220	Jordan	780	860	1040	500	480	490	600
FUL100250	London	790	840	1040	500	480	490	600
FUL100130	Bristol**	790	840	1040	480	480	490	600
FUL100310	Somerset	780	780	1060	480	490	490	600
FUL100150	Devon**	790	840	1040	480	480	490	600

**Alternative leg options available upon request.



Waterloo
Lounge Chair



Dublin
Highback Armchair



Rosie Cool
Highback Armchair



Drubin
Highback Armchair



Sutton
Highback Armchair



Penny
Highback Armchair



Padbury
Highback Armchair



Norfolk
Highback Armchair



Epping
Highback Armchair



Waterloo
2 Seater Lounge



Rosebud
2 / 2.5 Seater Lounge



Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100370	Waterloo	750	700	975	510	510	530	500
FUA500100	Dublin	660	790	1105	500	535	500	640
FUA500110	Rosie Cool	760	660	1110	500	530	510	635
FUA500130	Padbury*	600	700	1120	530	490	480	700
FUA500140	Norfolk	700	800	1110	490	480	490	640
FUA500120	Epping	610	800	1040	450	500	510	610

*Project only item: minimum order quantity required.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA500150	Drubin	760	670	1125	490	520	510	630
FUA500160	Sutton**	700	840	1050	510	490	490	600
FUA500170	Penny	770	660	1105	490	530	500	630
FUL100380	Waterloo - 2 Seater	1700	740	970	1495	510	520	500
FUL100360	Rosebud - 2 Seater	1410	750	960	1240	520	540	510
FUL100365	Rosebud - 2.5 Seater	1700	750	960	1530	520	540	510

*Project only item: minimum order quantity required. **Alternative leg options available upon request.



Leeton
2 Seater Lounge



Somerset
2 Seater Lounge



Sapphire
2 Seater Lounge



Jordan
2 Seater Lounge



Devon
2 Seater Lounge



Kew Mid Back
2 Seater Lounge



Kew Highback
2 Seater Lounge



London
2 Seater Lounge



Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100280	Leeton - 2 Seater	1200	780	800	1000	520	490	370
FUL100320	Somerset - 2 Seater	1300	780	1060	1000	490	490	600
FUL100160	Devon - 2 Seater**	1300	840	1040	1000	480	490	600
FUL100190	Kew Mid Back - 2 Seater	1200	840	940	1000	490	490	530

**Alternative leg options available upon request.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100340	Sapphire - 2 Seater	1200	740	880	1000	490	490	400
FUL100230	Jordan - 2 Seater	1250	860	1040	1000	480	490	600
FUL100210	Kew Highback - 2 Seater	1200	840	1040	1000	490	490	600
FUL100260	London - 2 Seater	1300	840	1040	1000	480	490	600



Cornwall
2 Seater Lounge



Bristol
2 Seater Lounge



Lawson
2 Seater Lounge



Lawson
3 Seater Lounge



Surrey
2 Seater Lounge



Wentworth
2 Seater Lounge



Wentworth
3 Seater Lounge



Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100110	Cornwall - 2 Seater**	1300	840	1040	1000	480	490	600
FUL100140	Bristol - 2 Seater**	1300	840	1040	1000	480	490	600
FUL100300	Surrey - 2 Seater	1200	840	1050	1000	490	490	600
FUL100390	Wentworth - 2 Seater	1600	840	840	1200	500	490	350

**Alternative leg options available upon request.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUL100410	Lawson - 2 Seater	1600	870	890	1200	500	490	400
FUL100420	Lawson - 3 Seater	2200	870	890	1800	500	490	400
FUL100400	Wentworth - 3 Seater	2200	840	840	1800	500	490	350



Armchairs

All Armchairs can be manufactured in customised upholstery and frame finishes. See pages 4 to 7.
 Images shown are for illustrative purposes only and may vary from custom options available.



Saki Armchair

Megan Armchair

Kelly Armchair

Adele Armchair



Bella Armchair

Willow Armchair

Mystique Armchair

Dane Armchair



Bev-a-Lee Armchair

Dorset Armchair

Urban Armchair

Mery Armchair

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100400	Bella	740	760	955	480	530	510	500
FUA100410	Willow	660	630	930	510	495	520	440
FUA100420	Mystique	660	630	915	510	480	525	400
FUA100440	Dane	585	550	920	480	440	470	470

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100450	Saki	580	590	810	465	460	470	360
FUA100460	Megan	710	550	735	460	440	475	300
FUA100470	Kelly	580	500	760	480	460	450	330
FUA100530	Adele*	610	600	850	450	440	470	420
FUA100480	Bev-a-Lee	660	640	760	480	450	490	360
FUA100490	Dorset	600	640	890	440	430	520	420
FUA100540	Urban	670	725	870	530	480	490	520
FUA100500	Mery*	725	600	770	490	460	460	370

*Project only item: minimum order quantity required.



Juliana Wide Armchair



Harmony Armchair



Diana Armchair



Stella Armchair

- Aluminium Frame -



Regina R Armchair



Regina G Armchair



Becca Armchair

- Aluminium Frame -



Jenary Armchair

- Aluminium Frame -



Coventry Armchair

- Aluminium Frame -



Gala Armchair

- Aluminium Frame -



Zenith Armchair

- Aluminium Frame -



Karol R Armchair



Karol G Armchair

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100390	Juliana - Wide	640	660	840	550	500	480	400
FUA100380	Harmony	620	640	1050	510	490	460	640
FUA100430	Diana	580	630	880	490	450	460	470
FUS100500	Becca	610	560	840	460	440	450	410
FUA100120	Jenary	580	600	950	500	450	470	510
FUA100100	Coventry	560	550	840	480	430	470	400
FUA100110	Gala	570	620	910	460	480	460	480

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100130	Stella	630	570	820	500	460	460	390
FUA100200	Regina R	560	550	860	480	430	470	430
FUA100210	Regina G	560	550	860	480	430	470	430
FUA100260	Zenith	590	540	950	480	480	470	420
FUA100220	Karol R	520	510	910	450	430	480	440
FUA100230	Karol G	520	510	910	450	430	480	440



Galaxy Armchair



Artella Armchair



Laura Armchair



Florence Armchair



Georgia Armchair



Danica Armchair



Vanessa Armchair



Valencia Armchair



Avenue Armchair



Froya Armchair



Geneva Armchair



Sienna Armchair

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100270	Galaxy - Black Steel frame	550	540	900	490	450	460	380
FUA100280	Galaxy - Chrome Steel Frame	550	540	900	490	450	460	380
FUA100290	Artella	565	570	890	485	445	500	430
FUA100300	Laura	600	580	920	490	475	500	450
FUA100190	Vanessa	570	530	890	490	440	480	420
FUA100320	Valencia	580	570	900	490	440	510	440
FUA100250	Avenue	580	600	870	490	490	470	450

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100140	Florence	600	540	860	480	470	460	440
FUA100240	Georgia	520	600	870	460	440	470	380
FUA100150	Dancia	550	570	910	470	470	460	460
FUA100160	Froya	590	530	830	480	440	470	420
FUA100170	Geneva	570	600	890	510	460	480	440
FUA100180	Sienna	560	590	860	480	480	490	360

Armchairs



Cairo Armchair

Paris Armchair

Delta Armchair

Tania Armchair



Manuela Armchair

Manuela Split Armchair

Joanne Armchair

Camille Armchair

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUA100330	Cairo	590	565	900	480	470	500	430
FUA100340	Paris	600	610	980	445	470	510	490
FUA100370	Delta	610	585	1010	490	460	500	545
FUA100310	Tania	570	575	810	490	470	510	320
FUA100350	Manuela	590	565	865	435	430	470	400
FUA100360	Manuela Split	590	585	865	435	450	470	400
FUA100520	Joanne*	620	530	760	500	440	450	320
FUA100510	Camille*	600	530	770	470	440	470	430

*Project only item: minimum order quantity required.



Side Chairs

All Side Chairs can be manufactured in customised upholstery and frame finishes. See pages 4 to 7. Images shown are for illustrative purposes only and may vary from custom options available.



Frankie Side Chair

Karol Side Chair

Regina Side Chair

Georgia Side Chair

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)
FUS100400	Frankie	525	530	925	445	400	510	430
FUS100100	Karol	460	510	910	450	430	480	440
FUS100200	Regina	470	560	850	460	440	460	440
FUS100300	Georgia	500	600	870	460	440	470	380



Day Chairs

Day Chairs are held in stock and available immediately in the finishes shown under each item. Custom Finishes are also available for some items where indicated.

KEY

CUSTOM FINISHES Customisable Upholstery



Classic Lowback Day Chair

- Slate Vinyl
- Champagne Vinyl
- Cream Frame



Classic Highback Day Chair

- Slate Vinyl
- Champagne Vinyl
- Cream Frame



Adjustable Highback Day Chair

- Ink Vinyl
- Latte Vinyl
- Cream Frame



Pressure Reducing Adjustable Highback Day Chair

- Ink Multistretch
- Cream Frame

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP208005 - Slate CHP208000 - Champagne	Classic Lowback	580	500	865 - 1015	460	440	450 - 605	440	160
CHP208055 - Slate CHP208050 - Champagne	Classic Highback	645	600	1160 - 1312	490	490	425 - 610	720	160
CHP208025 - Ink CHP208030 - Latte	Adjustable Highback	575 - 675	760 - 860	1070 - 1220	510 - 610	475 - 590	410 - 530	740 - 790	160
CHP208035 - Ink	Pressure Reducing Highback	575 - 675	760 - 860	1070 - 1220	510 - 610	475 - 590	410 - 530	740 - 790	160

Day Chairs



Maxi Adjustable Highback Day Chair

- Burgundy Vinyl
- Cream Frame



BC2 Waterfall Highback Day Chair

- Ink Vinyl
- Latte Vinyl
- Oatmeal
- Mocha
- Cream Frame
- CUSTOM FINISHES



BC4 Highback Day Chair

- Ink Vinyl
- Latte Vinyl
- Cream Frame



BC1 Lowback Day Chair

- Ink Vinyl
- Latte Vinyl
- Cream Frame
- CUSTOM FINISHES



BC1 Highback Day Chair

- Ink Vinyl
- Latte Vinyl
- Cream Frame
- CUSTOM FINISHES



BC1 Kingsize Highback Day Chair

- Ink Vinyl
- Cream Frame
- CUSTOM FINISHES

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP208010 - Burgundy	Maxi Adjustable Highback	800 - 1000	870 - 1000	1080 - 1470	600 - 800	420 - 510	440 - 580	750 - 1000	310
CHP196799 - Ink CHP196800 - Latte CHP208075 - Oatmeal CHP208076 - Mocha	BC2 Waterfall Highback	730	740	1020 - 1190	510	490	400 - 550	680	160
CHP196810 - Ink CHP196811 - Latte	BC4 Highback	630	780	990 - 1170	510	440 - 550	440 - 590	600	160
CHP196700 - Ink CHP196701 - Latte	BC1 Lowback	730	740	840 - 990	510	490	370 - 520	480	160
CHP196500 - Ink CHP196501 - Latte	BC1 Highback	730	740	990 - 1140	510	490	370 - 520	650	160
CHP196510 - Ink	BC1 Kingsize Highback	820	740	990 - 1140	600	490	370 - 520	650	200



Recliners

Recliners are held in stock and available immediately in the finish options shown for each item. Custom Finishes are also available for some items where indicated.

KEY

Customisable Upholstery

Motorised Product



Raphael Quattro Lift Recliner

Da Vinci Lift Recliner

Rembrandt Lift Recliner

- Bronze Microsuede Leather-look
- Marina Grey Microsuede Leather-look
- Nordic Stone Vinyl
-

- Gunmetal Microsuede
- Latte Microsuede

- Canyon Walnut Microsuede

Code	Description	Size	Overall Width (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP227290 - Bronze	Raphael Quattro	Small	740	1060	480	480	470	630	130
CHP227300 - Bronze CHP227301 - Marina Grey CHP227302 - Nordic Stone	Raphael Quattro	Medium	780	1090	520	500	470	680	130
CHP227310 - Bronze CHP227311 - Marina Grey CHP227312 - Nordic Stone	Raphael Quattro	Large	810	1080	550	530	500	700	130
CHP227319 - Bronze CHP227322 - Marina Grey	Raphael Quattro	Maxi	930	1050	650	540	510	650	180
CHP227510 - Gunmetal CHP227520 - Fog CHP227500 - Latte CHP227530 - Dove	Da Vinci	Medium	870	1050	530	500	480	630	159
CHP227550	Rembrandt	Maxi	1015	1095	560	535	495	600	159



Matisse Lift Recliner

- Santa Barbara Fog Microsuede



Renoir Lift Recliner

- Dusk Microsuede



Monet Lift Recliner

- Canyon Ocean Microsuede



Signature Space Saver Lift Recliner

- Heritage Vinyl
- Mink Microsuede
- Lagoon Microsuede



Signature Dual Action Lift Recliner

- Heritage Vinyl
- Mink Microsuede
- Lagoon Microsuede



Verona Space Saver Lift Recliner

- Mushroom Quilted Upholstery

Code	Description	Size	Overall Width (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP227560	Matisse	Large	965	1067	571	546	500	730	159
CHP227580	Renoir	Medium	800	1050	550	550	450	710	130
CHP227540	Monet	Large	890	1065	560	535	495	570	159
CHP197700 - Mink CHP197705 - Lagoon	Signature Space Saver	Small	980	740	480	460 - 515	480	508	145
CHP197710 - Mink CHP197715 - Lagoon CHP197716 - Heritage	Signature Space Saver	Medium	1050	780	510	490 - 515	480	620	145
CHP197720 - Mink CHP197725 - Lagoon CHP197730 - Heritage	Signature Space Saver	Large	1110	780	540	510 - 535	500	650	155
CHP197735 - Mink CHP197740 - Lagoon	Signature Dual Action	Small	980	740	480	460 - 515	480	508	145
CHP197745 - Mink CHP197750 - Lagoon CHP197751 - Heritage	Signature Dual Action	Medium	1050	780	510	490 - 515	480	620	145
CHP197755 - Mink CHP197760 - Lagoon CHP197765 - Heritage	Signature Dual Action	Large	1110	780	540	510 - 535	500	650	155
CHP227100 - Mushroom	Verona Space Saver	Medium	780	1100	540	500	500	630	130



Maui Dual Action Lift Recliner

Mocha Quilted Upholstery



Montana Single Action Maxi Lift Recliner

Bronze Microsuede Leather-look



Porto Zero Gravity Petite Lift Recliner

Fawn Microsuede



Oregon Lift Recliner

Fawn Microsuede



Idaho Lift Recliner

Fawn Microsuede



Altitude Vertical Lift Recliner

Bitter Chocolate Multistretch Beige Vinyl



Altitude Single Lift Recliner

Bitter Chocolate Multistretch Blue Vinyl



Terrigal Lift Recliner

Beige Vinyl Brown Vinyl



Nordic Mobicline Recliner

Ink Vinyl Latte Vinyl



Cove Recliner

Grey Multistretch Andreas Pumice Warwick
Navy Multistretch Bronx Steel Warwick



Prime Lift Recliner

Please enquire for fabric and vinyl options



PostureFit Recliner

Eastwood Ebony Multistretch Extreme Black Vinyl

Code	Description	Size	Overall Width (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP227110 - Mocha	Maui Dual Action	Medium	780	1100	540	500	500	630	130
CHP227400 - Bronze	Montana Single Action	Maxi	930	1050	650	540	510	650	220
CHP227200 - Fawn	Porto Zero Gravity	Petite	770	1120	490	500	490	690	150
CHP227150 - Charcoal CHP227155 - Latte	Oregon	Medium	800	1100	480	480	480	620	150
CHP227160 - Charcoal CHP227165 - Latte	Idaho	Medium	800	1100	480	480	480	620	150
CHP227700	Altitude Vertical	-	710	1200-1340	510	520	560 - 700	710	150

Code	Description	Size	Overall Width (mm)	Overall Height (mm)	Seat Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	SWL (kg)
CHP227705	Altitude Single	-	710	1200	510	520	560	710	150
CHP227900 - Beige CHP227910 - Dark Brown	Terrigal Electric	-	900	1120	500	490	490	630	140
CHP194715 - Ink CHP194700 - Latte	Nordic Mobicline	-	690	1120	530	490	530	680	140
CHP198960 - Grey / Andreas Pumice CHP198965 - Navy / Bronx Steel	Cove	Medium	775	1210	530	520 - 590	530	710	160
CHP198970 - Grey / Andreas Pumice CHP198975 - Navy / Bronx Steel	Cove	Large	825	1285	580	520 - 590	540	770	160
Range Available	Prime	-	710	1080	500	470	440	680	130
CHP198250 -	PostureFit	Medium	810 - 860	1190	510 - 560	420 - 480	490	800	185
CHP198260 -	PostureFit	Large	860 - 910	1210	560 - 610	420 - 480	490	800	185

Tables

- COFFEE TABLES
- SIDE TABLES
- DINING TABLES











Coffee Tables

All Coffee Tables may be customised to any of the premium or select stains available. See page 5.

Images shown are for illustrative purposes only and may vary from custom options available.

KEY

AVAILABLE CONFIGURATIONS

	Large Square		Large Round		Large Rectangle		Large Oval
	Square		Round		Rectangle		Oval



Euro Coffee Table



Peniche Coffee Table



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT100100	Euro Coffee Table - Oval	1195	590	500
FUT100110	Peniche Coffee Table - Round	500 ø	-	500

Coffee Tables



Piccolo Coffee Table Range



Delfi Coffee Table



Lily Coffee Table Range



Poli Coffee Table Range



Space Coffee Table Range



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT100140	Piccolo Coffee Table - Round	900 ø	-	450
FUT100130	Piccolo Coffee Table	1200	1200	330
FUT100170	Delfi Coffee Table	1200	1200	390
FUT100210	Lily Coffee Table	500	500	450
FUT100200	Lily Coffee Table	600	600	450
FUT100190	Lily Coffee Table - Round	500 ø	-	450
FUT100220	Lily Coffee Table	1200	600	450
FUT100160	Poli Coffee Table - Round	500 ø	-	465
FUT100150	Poli Coffee Table	500	500	455
FUT100230	Space Coffee Table - Round (Laminate)	600 ø	-	550
FUT100250	Space Coffee Table - Round (Veneer)	600 ø	-	550
FUT100240	Space Coffee Table - Round (Laminate)	1100 ø	-	550
FUT100260	Space Coffee Table - Round (Veneer)	1100 ø	-	550



Side Tables

All Side Tables may be customised to any of the premium or select stains available. See page 5.

Images shown are for illustrative purposes only and may vary from custom options available.

KEY

AVAILABLE CONFIGURATIONS

Large Round (grey circle) Round (small grey circle)



Piccolo Side Table



Space Side Table



Euro Side Table



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT300110	Piccolo Side Table	450 ø	-	600
FUT300120	Space Side Table (Laminate)	450 ø	-	550
FUT300130	Space Side Table (Veneer)	450 ø	-	550
FUT300140	Euro Side Table	595 ø	-	600

Side Tables



Alta Side Table



Peniche Side Table



Lily Side Table



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT300100	Alta Side Table	400 ø	-	600
FUT300170	Peniche Side Table	500 ø	-	600
FUT300150	Lily Side Table	500 ø	-	550




Dining Tables

All Dining Tables may be customised to any of the premium or select stains available. See page 5.

Images shown are for illustrative purposes only and may vary from custom options available.

KEY

AVAILABLE CONFIGURATIONS



Dallas Dining Table



Faison Dining Table



Solo Dining Table



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT500180	Dallas Dining Table	1200 ø	-	760
FUT500350	Faison Dining Table	1200 ø	-	760
FUT020260	Solo Dining Table	900 ø	-	750



Lily Dining Table Range



Tetis Dining Table Range



Frankie Dining Table



City Dining Table Range



Signature Dining Table Range



Peniche Dining Table Range



Via Dining Table



Clover Leaf Dining Table



Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT500200	Lily Dining Table	1000	1000	760
FUT500220	Lily Dining Table	1600	1000	760
FUT500250	Tetis Dining Table	900	900	760
FUT500260	Tetis Dining Table	1500	900	760
FUT500240	Tetis Dining Table - Round	1000 ϕ	-	760
FUT500140	Signature Dining Table	1000	1000	750
FUT500150	Signature Dining Table (Laminate)	1000	1000	750
FUT500120	Signature Dining Table	1800	900	750
FUT500130	Signature Dining Table (Laminate)	1800	900	750
FUT500160	Signature Dining Table - Round	1000 ϕ	-	750
FUT500170	Signature Dining Table - Round (Laminate)	1000 ϕ	-	750
FUT500330	Peniche Dining Table	1000	1000	760
FUT500310	Peniche Dining Table	1600	1000	760

Code	Description	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)
FUT500100	Frankie Dining Table - Oval	2000	1000	760
FUT500360	City Dining Table	800	800	750
FUT500370	City Dining Table	900	900	750
FUT500380	City Dining Table	1000	1000	750
FUT500390	City Dining Table	1200	800	750
FUT500400	City Dining Table	1500	800	750
FUT500410	City Dining Table	1800	800	750
FUT500420	City Dining Table	1800	900	750
FUT500430	City Dining Table	2400	900	750
FUT500190	Via Dining Table	900	900	750
FUT500110	Clover Leaf Dining Table	1200	1200	760

Storage

- SIDEBOARDS & BUFFETS
- CABINETS
 - STORAGE CABINETS
 - BEDSIDE CABINETS
 - TV AND MEDIA CABINETS



Sideboards & Buffets

All Sideboards and Buffets may be customised to any of the premium or select stains available. See page 5.

Images shown are for illustrative purposes only and may vary from custom options available.



Camden Sideboard



Alice Sideboard



Allison Sideboard



Shadie Sideboard



Lily Sideboard

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
FUB100130	Camden Sideboard - 2 Module	1200	500	800
FUB100110	Alice Sideboard - 3 Module	1500	450	900
FUB100120	Allison Sideboard	1200	400	760
FUB100100	Shadie Sideboard	1200	450	750
FUB100140	Lily Sideboard	1400	500	760



Cabinets

All Cabinets may be customised to any of the premium or select stains available. See page 5.

Images shown are for illustrative purposes only and may vary from custom options available.



Belmont Bedside Cabinet



Freya Bedside Cabinet



Freya Premium Bedside Cabinet



Hansel Bedside Cabinet



Duke Mini Bedside Cabinet



Duke Wide Bedside Cabinet



Acute Bedside Cabinet



Regent Bedside Cabinet



With Pull Out Top Shelf

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
BEA826230	Belmont Bedside Cabinet (Veneer)	470	450	715
BEA826360	Belmont Bedside Cabinet (Laminate)	470	450	715
BEA826260	Freya Bedside Cabinet	495	450	700
BEA826280	Freya Premium Bedside Cabinet (Veneer)	500	450	700
BEA826270	Freya Premium Bedside Cabinet (Laminate)	500	450	700

All Bedside Cabinets come standard with built in lock.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
BEA826290	Hansel Bedside Cabinet	495	450	700
BEA826210	Duke Mini Bedside Cabinet (Veneer)	500	450	700
BEA826200	Duke Mini Bedside Cabinet (Laminate)	500	450	700
BEA826220	Duke Wide Bedside Cabinet (Veneer)	570	450	750
BEA826380	Acute Bedside Cabinet	465	525	840
BEA826250	Regent Bedside Cabinet (Veneer)	475	450	705

All Bedside Cabinets come standard with built in lock.



With Hutch



With Hutch

Austin Bedside Cabinet

Belinda Bedside Cabinet




Alice Storage Cabinet

Camden Storage Case


Austin Wardrobe

Belinda Wardrobe

- Available in Slim -



With Hutch



With Hutch

Chelsea Bedside Cabinet

Derik Bedside Cabinet



Camden Media Cabinet

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
BEA826300	Austin Bedside Cabinet	465	450	810
BEA826305	Austin Bedside Cabinet + Hutch	465	450	1020
BEA826310	Belinda Bedside Cabinet	465	450	810
BEA826315	Belinda Bedside Cabinet + Hutch	465	450	1020
BEA826320	Chelsea Bedside Cabinet	465	450	810
BEA826330	Chelsea Bedside Cabinet + Hutch	465	450	1020
BEA826340	Derik Bedside Cabinet	465	450	810
BEA826350	Derik Bedside Cabinet + Hutch	465	450	1020

All Bedside Cabinets come standard with built in lock.

Code	Description	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
FUL053110	Alice Storage Cabinet	885	400	1900
BEA826510	Camden Storage Case	900	350	1800
BEA826520	Austin Slim Wardrobe	465	570	1800
BEA826530	Austin Wardrobe	605	570	1800
BEA826540	Belinda Wardrobe	655	570	1800
FUM100100	Camden Media Cabinet	1800	500	500

All Bedside Cabinets come standard with built in lock.

Beds

aidacare
FL250
Low height, high style

With its impressive low height of 110mm, injury to the resident can be significantly reduced or avoided in the event of a fall from bed, therefore reducing the need for side-rails.



- TO SCAN QR CODE
- ① Open phone camera
 - ② Hover over QR code
 - ③ Click on link at the top



Additional Head and Footboard range available



Tallo



Eyre



Como



Sophia


aspire
ComfiMotion®
Collection

We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion* Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.



- TO SCAN QR CODE
- ① Open phone camera
 - ② Hover over QR code
 - ③ Click on link at the top




CARE

LUXE



◀ Companion Bed

Headboards



Chesterfield



Allure

Mattress Options



Pocket Spring



Memory Foam

Alfresco Furniture



ALFRESCO Furniture



WARRANTY

Subject to any express or implied statutory requirements, Aidacare provides the following warranty on custom furniture:

5 YEAR WARRANTY
INTERNAL FURNITURE

3 YEAR WARRANTY
EXTERNAL FURNITURE

- Structural Warranty refers to furniture frames and excludes moving hardware components and fixings which have a 1 Year Warranty.
- Warranty is offered on the basis that the warranted product is used for the purpose of which it was designed.
- Warranty covers defects in material and workmanship - general wear-and-tear and/or misuse, natural variations in veneers, timber and laminate are not considered defects.
- Warranty does not apply to damage cause by accident, neglect and improper operation.



COPYRIGHT

© Copyright Aidacare Pty Ltd 2021 Version 1.3

No part of this catalogue may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means.

DISCLAIMER

Due care has been taken to ensure all the information contained in this catalogue is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. We reserve the right to modify the design or appearance of any product displayed in this catalogue without prior notice.



PROJECTS BY AIDACARE DISPLAY CENTRE

Consult • Design • Project Manage

The Projects by Aidacare Display Centre is a dedicated, visually-tailored solution for Furniture, Fittings & Equipment. A one-stop solution for Aged Care Greenfield & Brownfield developments and Commercial projects.



projectsbyaidacare.com.au

To see more, visit our website
atamaliving.com.au



TO SCAN QR CODE

- ① Open phone camera
- ② Hover over QR code
- ③ Click on link at the top



Solutions for

- Bedroom Suites
- Consultation Rooms
- Sluice Rooms
- Ceiling Hoist Solutions
- Interior Design
- Common & Living Spaces
- Furniture, Fittings & Equipment





ACA44-23 060622

P 1300 083 661
E info@atamaliving.com.au
W atamaliving.com.au

ATAMA  LIVING
by aidacare